


Anti-Spam1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P38567

Item Number: CM0023218



Introduction

Rabbit polyclonal antibody targeting human SPAM1 (PH20) hyaluronidase. Validated for WB, IF-P, IHC-P, ELISA with cross-reactivity to mouse and rat. Suitable for reproductive biology and fertility research.

[Learn More](#)

Parameter	Value
Product Name	SPAM1 Rabbit pAb
Remarks/Alias	HYA1, PH20, HYAL1, HYAL3, HYAL5, PH-20, SPAG15, HEL-S-96n, SPAM1
Species	Human
Gene ID	6677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human SPAM1 (NP_003108.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38567
Protein Weight	58 kDa
Shipping	Ice bag