

Anti-Sparc Rabbit Polyclonal Antibody For Ihc-P, Elisa - P09486

Item Number: CM0020698



Introduction

High-quality rabbit polyclonal antibody against human SPARC (Osteonectin, BM-40). Validated for immunohistochemistry (IHC-P) and ELISA; suitable for tissue staining and quantitative assays in human, mouse, and rat samples. Ideal for extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	SPARC Rabbit pAb
Remarks/Alias	ON; ONT; OI17; BM-40; SPARC
Species	Human
GeneID (Human)	6678
GeneID	6678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human SPARC (NP_003109.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09486
Protein Weight	35kDa
Shipping	Ice bag