


Anti-Sparcl1 Rabbit Polyclonal Antibody For Wb And Elisa - Q14515

Item Number: CM0021840



Introduction

High-quality polyclonal antibody against human SPARCL1 (Hevin, MAST9). Validated for WB and ELISA, reactive with human, mouse, and rat. Suitable for extracellular matrix and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SPARCL1 Rabbit pAb
Remarks/Alias	SC1; MAST9; PIG33; MAST 9; SPARCL1
Species	Human
GeneID (Human)	8404
GeneID	8404
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-300 of human SPARCL1 (NP_004675.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14515
Protein Weight	75kDa
Shipping	Ice bag