

Anti-Spata4 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ney3

Item Number: CM0029587

Introduction

SPATA4 Rabbit Polyclonal Antibody validated for WB and ELISA, recognizing human, mouse, and rat SPATA4 (35kDa). Suitable for apoptosis and spermatogenesis research in testis tissue.

[Learn More](#)

Parameter	Value
Product Name	SPATA4 Rabbit pAb
Remarks/Alias	FAP178; SPEF1B; TSARG2; CFAP178; SPATA4
Species	Human
GeneID	132851
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human SPATA4 (NP_653245.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NEY3
Protein Weight	35kDa
Shipping	Ice bag