


Anti-Spata5 Polyclonal Antibody For Wb, Elisa - Q8Nb90

Item Number: CM0022623



Introduction

Rabbit polyclonal antibody targeting human SPATA5 (AFG2, SPAF), validated for WB and ELISA. Detects human, mouse, and rat SPATA5; relevant for studies on chromatin replication, ribosome biogenesis, and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	SPATA5 Rabbit pAb
Remarks/Alias	AFG2; SPAF; EHLMRS; NEDHSB; SPATA5
Species	Human
GeneID (Human)	166378
GeneID	166378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 228-304 of human SPATA5 (NP_001304728)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB90
Protein Weight	98kDa
Shipping	Ice bag