

Anti-Spc24 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nbt2

Item Number: CM0024675



Introduction

Validate SPC24's role in kinetochore assembly and chromosome segregation with this rabbit polyclonal antibody. Suitable for WB and ELISA; cross-reacts with human and mouse SPC24.

[Learn More](#)

Parameter	Specification
Product Name	SPC24 Rabbit pAb
Remarks/Alias	SPBC24; SPC24
Species	Human
GeneID (human)	147841
GeneID	147841
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human SPC24 (NP_872319.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NBT2
Protein Weight	22kDa
Shipping	Ice bag