

Anti-Spc25 Polyclonal Antibody For Wb, Elisa - Q9Hbm1

Item Number: CM0015953



Introduction

Rabbit polyclonal antibody targeting human SPC25 (kinetochore protein Spc25), validated for Western blot and ELISA applications. Suitable for research on chromosome segregation, mitosis, and kinetochore function. Swiss-Prot Q9HBM1.

[Learn More](#)

Parameter	Value
Product Name	SPC25 Rabbit pAb
Remarks/Alias	AD024; SPBC25; hSpc25; SPC25
Species	Human
GeneID (Human)	57405
GeneID	57405
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-224 of human SPC25 (NP_065726.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HBM1
Protein Weight	26kDa
Shipping	Ice bag