

Anti-Spdef Polyclonal Antibody For Wb, Elisa - 095238

Item Number: CM0023936



Introduction

Rabbit polyclonal antibody against human SPDEF, validated by WB and ELISA. Detects a 38 kDa band and is suitable for prostate cancer research and prostate-specific antigen promoter studies.

[Learn More](#)

Parameter	Value
Product Name	SPDEF Rabbit pAb
Remarks/Alias	PDEF; bA375E1.3; SPDEF
Species	Human
Gene ID (Human)	25803
Gene ID	25803
Immunogen	Synthetic peptide, aa 1-100 of human SPDEF (NP_036523.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95238
Protein Weight	38kDa
Shipping	Ice bag