

Anti-Specc1L Polyclonal Antibody For Wb, Ihc-P, Elisa - Q69Yq0

Item Number: CM0019478

Introduction

High-quality rabbit polyclonal antibody against human SPECC1L (Cytospin-A). Validated for WB, IHC-P, ELISA. Ideal for studying cytokinesis, spindle organization, and cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	SPECC1L Rabbit pAb
Remarks/Alias	TBHS; CYTSA; GBBB2; TBHS1; OBLFC1; SPECC1L
Species	Human
Gene ID	23384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-500 of human SPECC1L (NP_056145.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q69YQ0
Protein Weight	125kDa
Shipping	Ice bag