

Anti-Spem1 Polyclonal Antibody For Wb, Elisa - Q8N4L4

Item Number: CM0028209

Introduction

High-quality rabbit polyclonal antibody targeting human SPEM1 (C17orf83), validated for WB and ELISA. Detects a membrane/cytoplasmic protein essential for spermatid cytoplasm removal during spermatogenesis. Cross-reacts with rat.

[Learn More](#)

Parameter	Value
Product Name	SPEM1 Rabbit pAb
Remarks/Alias	C17orf83; SPEM1
Species	Human
Gene ID (Human)	374768
Gene ID	374768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 203-309 of human SPEM1 (NP_955371.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8N4L4
Protein Weight	35kDa
Shipping	Ice bag