


Anti-Spermine Synthase Monoclonal Antibody For Wb, Elisa - P52788

Item Number: CM0011767



Introduction

Rabbit monoclonal antibody targeting human Spermine synthase (SMS), suitable for WB and ELISA applications. Detects human, mouse, and rat species. Ideal for studying polyamine biosynthesis and related disorders.

[Learn More](#)

Parameter	Value
Product Name	Spermine synthase Rabbit mAb
Remarks/Alias	SRS; SpS; MRSR; SPMSY; MRXSSR; Spermine synthase
Species	Human
Gene ID	6611
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Spermine synthase (P52788).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52788
Protein Weight	41 kDa
Shipping	Ice bag