


Anti-Spermine Synthase (Sms) Rabbit Polyclonal Antibody For Wb And Elisa - P52788

Item Number: CM0027153



Introduction

High-quality rabbit polyclonal antibody targeting human spermine synthase (SMS). Validated for Western blot and ELISA applications. Cross-reacts with human and mouse. Essential for polyamine biosynthesis and neurodevelopment research.

[Learn More](#)

Parameter	Value
Product Name	SMS Rabbit pAb
Remarks/Alias	SRS, SpS, MRSR, SPMSY, MRXSSR, SMS
Species	Human
GeneID (Human)	6611
GeneID	6611
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-366 of human SMS (NP_004586.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52788
Protein Weight	41kDa
Shipping	Ice bag