


Anti-Sphk1 Polyclonal Antibody For Ihc-P, Elisa - Q9Nya1

Item Number: CM0020191



Introduction

Rabbit polyclonal antibody to human sphingosine kinase 1 (SPHK1) for IHC-P and ELISA, cross-reactive with mouse and rat. Detects a 43kDa band. Ideal for researching S1P-mediated signaling, inflammation, and neuroinflammation pathways.

[Learn More](#)

Parameter	Value
Product Name	SPHK1 Rabbit pAb
Remarks/Alias	SPHK
Species	Human
GeneID (Human)	8877
GeneID	8877
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 201-300 of human SPHK1 (NP_001136074.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYA1
Protein Weight	43kDa
Shipping	Ice bag