


## Anti-Sphk1 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Nya1

Item Number: CM0006796



### Introduction

SPHK1 rabbit monoclonal antibody validated for Western blot and ELISA. Targets human sphingosine kinase 1 (SWISS Q9NYA1), a lipid kinase involved in sphingolipid metabolism and inflammation.

[Learn More](#)

| Parameter        | Value                                                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SPHK1 Rabbit mAb                                                                                                                           |
| Remarks/Alias    | SPHK; SPHK1                                                                                                                                |
| Species          | Human                                                                                                                                      |
| GeneID (Human)   | 8877                                                                                                                                       |
| GeneID           | 8877                                                                                                                                       |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human SPHK1 (NP_001136074.1). |
| Source           | Rabbit                                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                                      |
| Application      | WB, ELISA                                                                                                                                  |
| Cross Reactivity | Human                                                                                                                                      |
| SWISS            | Q9NYA1                                                                                                                                     |
| Protein Weight   | 43kDa                                                                                                                                      |
| Shipping         | Ice bag                                                                                                                                    |