


Anti-Sphk1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nya1

Item Number: CM0018660



Introduction

Rabbit polyclonal antibody targeting human sphingosine kinase 1 (SPHK1), validated for Western blot (WB) and ELISA. Cross-reacts with mouse samples. Supports research in cell growth, inflammation, and lipid signaling.

[Learn More](#)

Parameter	Value
Product Name	SPHK1 Rabbit pAb
Remarks/Alias	SPHK; SPHK1
Species	Human
GeneID (Human)	8877
GeneID	8877
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-384 of human SPHK1 (NP_001136074.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NYA1
Protein Weight	43 kDa
Shipping	Ice bag