

Anti-Sphk2 Polyclonal Antibody For If/Icc And Elisa - Q9Nra0

Item Number: CM0028374

Introduction

SPHK2 Rabbit pAb is a polyclonal antibody against human Sphingosine kinase 2, validated for IF/ICC and ELISA. Ideal for investigating sphingolipid signaling, apoptosis, and mitochondrial function in human samples.

[Learn More](#)

Parameter	Value
Product Name	SPHK2 Rabbit pAb
Remarks/Alias	SK 2; SK-2; SPK 2; SPK-2
Species	Human
Gene ID (Human)	56848
Gene ID	56848
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-654 of human SPHK2 (NP_064511.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9NRA0
Protein Weight	69kDa
Shipping	Ice bag