

Anti-Sphk2 Polyclonal Antibody For Wb And Elisa - Q9Nra0

Item Number: CM0029742

Introduction

SPHK2 rabbit polyclonal antibody for WB and ELISA, targeting human sphingosine kinase 2 (SPHK2). Cross-reacts with human and mouse. Ideal for studying S1P signaling, apoptosis, mitochondrial respiration, and epigenetic gene regulation.

[Learn More](#)

Parameter	Value
Product Name	SPHK2 Rabbit pAb
Remarks/Alias	SK 2; SK-2; SPK 2; SPK-2; SPHK2
Species	Human
GeneID	56848
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human SPHK2 (NP_001191088.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NRA0
Protein Weight	69kDa
Shipping	Ice bag