

Anti-Sphk2 Polyclonal Antibody For Wb, Elisa - Q9Nra0

Item Number: CM0025108

Introduction

Rabbit polyclonal antibody targeting human SPHK2, validated for Western blot (WB) and ELISA applications. Detects human, mouse, and rat SPHK2 with a predicted molecular weight of 69 kDa. Ideal for apoptosis and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SPHK2 Rabbit pAb
Remarks/Alias	SK 2; SK-2; SPK 2; SPK-2; SPHK2
Species	Human
Gene ID (Human)	56848
Gene ID	56848
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human SPHK2 (NP_001191088.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRA0
Protein Weight	69kDa
Shipping	Ice bag