

Anti-Spice1 Polyclonal Antibody For Western Blot And Elisa - Q8N0Z3

Item Number: CM0029869

Introduction

Rabbit polyclonal antibody targeting human SPICE1 (Spindle and centriole-associated protein 1), validated for Western blotting and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studies on centriole duplication, spindle formation, and mitosis.

[Learn More](#)

Parameter	Description
Product Name	SPICE1 Rabbit pAb
Remarks/Alias	SPICE; CCDC52; SPICE1
Species	Human
GeneID (Human)	152185
GeneID	152185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SPICE1 (NP_653319.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N0Z3
Protein Weight	96kDa
Shipping	Ice bag