

Anti-Spike Glycoprotein Monoclonal Antibody For Wb, If/Icc, Elisa - P0Dtc2

Item Number: CM0002102



Introduction

Recombinant rabbit monoclonal antibody targeting SARS-CoV-2 spike glycoprotein (S protein) for WB, IF/ICC, ELISA applications. Detects a band of ~141 kDa. Suitable for viral entry and host interaction studies.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 3CLpro Rabbit mAb
Remarks/Alias	Spike glycoprotein; S; SARS-CoV-2 3CLpro
Species	Severe acute respiratory syndrome coronavirus 2
GeneID	43740578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of coronavirus 3CLpro (YP_009725295.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTC2
Protein Weight	141kDa
Shipping	Ice bag