


Anti-Spin1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y657

Item Number: CM0027696



Introduction

Rabbit polyclonal antibody targeting SPIN1 (Spindlin-1), a chromatin reader that recognizes H3K4me3K9me3 and activates Wnt signaling. Validated for WB, IHC-P, ELISA, with cross-reactivity to human, mouse, and rat. Ideal for cancer and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	SPIN1 Rabbit pAb
Remarks/Alias	SPIN; TDRD24; SPIN1
Species	Human
Gene ID (Human)	10927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-262 of human SPIN1 (NP_006708.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y657
Protein Weight	30kDa
Shipping	Ice bag