

## Anti-Spin2B Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bpz2

Item Number: CM0025428



### Introduction

Anti-SPIN2B Rabbit Polyclonal Antibody (29kDa) validated for WB and ELISA. Detects human SPIN2B and cross-reacts with mouse and rat. SPIN2B is a nuclear protein involved in cell cycle regulation and apoptosis, with H3K4me3-binding activity.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SPIN2B Rabbit pAb                                                                                                    |
| Remarks/Alias    | SPIN-2; TDRD26; SPIN-2B; dj323P24.2; SPIN2_duplicate; SPIN2B                                                         |
| Species          | Human                                                                                                                |
| Gene ID          | 474343                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human SPIN2B (NP_001006682.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | Q9BPZ2                                                                                                               |
| Protein Weight   | 29kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |