

Anti-Spink1 Rabbit Polyclonal Antibody For If/Icc, Elisa - P00995

Item Number: CM0014738



Introduction

Anti-SPINK1 rabbit polyclonal antibody validated for IF/ICC and ELISA applications. Detects human SPINK1, a serine protease inhibitor involved in pancreatic function and reproductive biology. Ideal for IVD and research use.

[Learn More](#)

Parameter	Value
Product Name	SPINK1 Rabbit pAb
Remarks/Alias	TCP, PCTT, PSTI, TATI, Spink3, SPINK1
Species	Human
GeneID (Human)	6690
GeneID	6690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-79 of human SPINK1 (NP_003113.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P00995
Protein Weight	9kDa
Shipping	Ice bag