

Anti-Spink5/Lekti Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nq38

Item Number: CM0002918



Introduction

High-quality Anti-SPINK5/LEKTI Rabbit Monoclonal Antibody validated for WB and ELISA. Cross-reacts with Human, Mouse, Rat. Recombinant immunogen, suitable for research on serine protease inhibitor and skin barrier function.

[Learn More](#)

Parameter	Value
Product Name	SPINK5 / LEKTI Rabbit mAb
Remarks/Alias	NS; NETS; LEKTI; LETKI; VAKTI; SPINK5 / LEKTI
Species	Human
GeneID (Human)	11005
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human SPINK5/LEKTI (Q9NQ38)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ38
Protein Weight	121kDa
Shipping	Ice bag