


Anti-Spint1 (Hai-1) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43278

Item Number: CM0027779



Introduction

Rabbit polyclonal antibody against human SPINT1 (HAI-1), a Kunitz-type protease inhibitor, validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat; recombinant immunogen; predicted weight 58 kDa.

[Learn More](#)

Parameter	Value
Product Name	SPINT1 Rabbit pAb
Remarks/Alias	HAI; HAI1; MANSC2; SPINT1
Species	Human
GeneID (Human)	6692
GeneID	6692
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 234-513 of human SPINT1 (NP_003701.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43278
Protein Weight	58kDa
Shipping	Ice bag