

## Anti-Spk (Topk) Rabbit Polyclonal Antibody - Q96Kb5

Item Number: CM0030411



### Introduction

Anti-SPK (TOPK) rabbit polyclonal antibody for WB, IHC-P, and ELISA. Reacts with human and rat samples. Detects a 36 kDa protein; ideal for cancer/testis antigen and kinase signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | SPK Rabbit pAb                                                                                                  |
| Remarks/Alias    | SPK; CT84; TOPK; HEL164; Nori-3                                                                                 |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 55872                                                                                                           |
| GeneID           | 55872                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-322 of human SPK (NP_060962.2). |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Rat                                                                                                      |
| SWISS            | Q96KB5                                                                                                          |
| Protein Weight   | 36kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |