

Anti-Spock3 Polyclonal Antibody For Wb, Elisa - Q9Bq16

Item Number: CM0018160

Introduction

Rabbit polyclonal antibody against human SPOCK3 (Testican-3) validated for WB, ELISA. Cross-reacts with human, mouse, rat. For neurogenesis and tumor invasion research.

[Learn More](#)

Parameter	Value
Product Name	SPOCK3 Rabbit pAb
Remarks/Alias	TES-3; TICN3; HSAJ1454; SPOCK3
Species	Human
Gene ID	50859
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human SPOCK3 (NP_001035249.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ16
Protein Weight	49 kDa
Shipping	Ice bag