


Anti-Spon2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bud6

Item Number: CM0014296



Introduction

High-quality rabbit polyclonal antibody targeting human SPON2 (MINDIN/DIL-1), validated in WB, IHC-P, and ELISA. Recognizes a secreted extracellular matrix protein with roles in cell adhesion and innate immunity. Cross-reacts with rat.

[Learn More](#)

Parameter	Specification
Product Name	SPON2 Rabbit pAb
Remarks/Alias	DIL1; DIL-1; MINDIN; M-SPONDIN; SPON2
Species	Human
Gene ID (Human)	10417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-150 of human SPON2 (NP_036577.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BUD6
Protein Weight	36kDa
Shipping	Ice bag