


Anti-Sppl2A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tct8

Item Number: CM0016126



Introduction

High-quality rabbit polyclonal antibody against human SPPL2A (IMP3/PSL2), validated for WB, IHC-P, ELISA. Detects endogenous SPPL2A in human, mouse, rat. Ideal for immunology and protease research.

[Learn More](#)

Parameter	Value
Product Name	SPPL2A Rabbit pAb
Remarks/Alias	IMP3, PSL2, IMD86, SPPL2A
Species	Human
Gene ID (Human)	84888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 421-520 of human SPPL2A (NP_116191.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TCT8
Protein Weight	58kDa
Shipping	Ice bag