


Anti-Spr (Sepiapterin Reductase) Polyclonal Antibody For Wb, Elisa - P35270

Item Number: CM0017478



Introduction

Rabbit polyclonal antibody against human SPR (Sepiapterin reductase), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying tetrahydrobiopterin biosynthesis and related dystonia research.

[Learn More](#)

Parameter	Value
Product Name	SPR Rabbit pAb
Remarks/Alias	SDR38C1; SPR
Species	Human
Gene ID (Human)	6697
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human SPR (NP_003115.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35270
Protein Weight	28kDa
Shipping	Ice bag