

Anti-Sprouty 2 Monoclonal Antibody For Wb, Ihc-P, Elisa - O43597

Item Number: CM0010433



Introduction

Rabbit monoclonal antibody targeting human Sprouty 2 (SPRY2) for Western blot, IHC-P, and ELISA applications, with cross-reactivity to rat. Suitable for studying FGF signaling and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	Sprouty 2 Rabbit mAb
Remarks/Alias	IGAN3; hSPRY2; Sprouty 2
Species	Human
GeneID (Human)	10253
GeneID	10253
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 216-315 of human Sprouty 2 (O43597).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O43597
Protein Weight	35kDa
Shipping	Ice bag