

Anti-Sprrr1A Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P35321

Item Number: CM0019595

Introduction

SPRR1A Rabbit Polyclonal Antibody (pAb) targets Cornifin-A (SPRR1A). Validated for WB, IF-P, ELISA with human, mouse, rat cross-reactivity. Ideal for keratinocyte differentiation, skin barrier, and cornified envelope research.

[Learn More](#)

Parameter	Value
Product Name	SPRR1A Rabbit pAb
Remarks/Alias	SPRK; SPRR1A
Species	Human
Gene ID (Human)	6698
Gene ID	6698
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-89 of human SPRR1A (NP_005978.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35321
Protein Weight	10kDa
Shipping	Ice bag