

Anti-Sprp3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubc9

Item Number: CM0015856

Introduction

Rabbit polyclonal antibody against human SPRR3, validated for Western blot, immunofluorescence/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for keratinocyte and cornified envelope research.

[Learn More](#)

Parameter	Specification
Product Name	SPRR3 Rabbit pAb
Remarks/Alias	SPRR3
Species	Human
GeneID (Human)	6707
GeneID	6707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human SPRR3 (NP_005407.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBC9
Protein Weight	17kDa
Shipping	Ice bag