


Anti-Sprtn Rabbit Polyclonal Antibody For Western Blot And Elisa - Q9H040

Item Number: CM0016549



Introduction

Rabbit polyclonal antibody targeting human SPRTN (DVC1/Spartan), a DNA-dependent metalloprotease critical for DNA-protein crosslink repair and genomic stability. Validated for Western blot and ELISA; detects a 55 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	SPRTN Rabbit pAb
Remarks/Alias	DVC1; PRO4323; spartan; C1orf124; SPRTN
Species	Human
Gene ID	83932
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human SPRTN (NP_001010984.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H040
Protein Weight	55kDa
Shipping	Ice bag