

Anti-Spry2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43597

Item Number: CM0023422



Introduction

SPRY2 Rabbit Polyclonal Antibody validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Reacts with human and mouse. Suitable for studying FGF signaling, development, and EGFR ubiquitination.

[Learn More](#)

Parameter	Value
Product Name	SPRY2 Rabbit pAb
Remarks/Alias	IGAN3; hSPRY2; SPRY2
Species	Human
GeneID (Human)	10253
GeneID	10253
Immunogen	Synthetic peptide [A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SPRY2 (NP_005833.1)]
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43597
Protein Weight	35kDa
Shipping	Ice bag