

Anti-Spry3 Rabbit Polyclonal Antibody For Wb, Elisa - 043610

Item Number: CM0016484

Introduction

Rabbit polyclonal antibody against human SPRY3, validated for WB and ELISA. Detects Protein sprouty homolog 3, a MAPK signaling antagonist involved in neurogenesis and organogenesis. Suitable for brain and heart tissue studies.

[Learn More](#)

Parameter	Value
Product Name	SPRY3 Rabbit pAb
Remarks/Alias	spry-3; SPRY3
Species	Human
Gene ID (Human)	10251
Gene ID	10251
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human SPRY3 (NP_005831.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43610
Protein Weight	31kDa
Shipping	Ice bag