


Anti-Spsb1 Polyclonal Antibody For Wb And Elisa - Q96Bd6

Item Number: CM0015199



Introduction

Rabbit polyclonal antibody against human SPSB1, suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Targets a key regulator of nitric oxide production and ubiquitin-mediated degradation.

[Learn More](#)

Parameter	Value
Product Name	SPSB1 Rabbit pAb
Remarks/Alias	SSB1; SSB-1; SPSB1
Species	Human
GeneID (human)	80176
GeneID	80176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human SPSB1 (NP_079382.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96BD6
Protein Weight	31kDa
Shipping	Ice bag