

Anti-Sptbn1 Rabbit Polyclonal Antibody For Wb And Elisa - Q01082

Item Number: CM0026569



Introduction

Rabbit polyclonal antibody against human SPTBN1 (beta-II spectrin), validated for WB and ELISA. Detects human and mouse SPTBN1, ~275kDa. Ideal for cytoskeleton and membrane research.

[Learn More](#)

Parameter	Value
Product Name	SPTBN1 Rabbit pAb
Remarks/Alias	ELF; SPTB2; DDISBA; HEL102; betaSplI; SPTBN1
Species	Human
Gene ID (Human)	6711
Gene ID	6711
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2100-2200 of human SPTBN1 (NP_003119.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q01082
Protein Weight	275kDa
Shipping	Ice bag