


Anti-Sptlc1 Polyclonal Antibody For Wb, Elisa - 015269

Item Number: CM0028890



Introduction

SPTLC1 Rabbit Polyclonal Antibody validated for Western blotting (WB) and ELISA. Detects human, mouse, rat SPTLC1. Ideal for sphingolipid biosynthesis and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SPTLC1 Rabbit pAb
Remarks/Alias	HSN1; LBC1; LCB1; SPT1; SPTI; ALS27; HSN1; SPTLC1
Species	Human
GeneID	10558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-310 of human SPTLC1 (NP_006406.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15269
Protein Weight	53kDa
Shipping	Ice bag