

## Anti-Sqle Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14534

Item Number: CM0020603



### Introduction

High-quality rabbit polyclonal antibody against human squalene monooxygenase (SQLE) for WB, IHC-P, IF/ICC, ELISA. Validated in human, mouse, rat. Ideal for cholesterol metabolism and liver disease research.

[Learn More](#)

| Parameter        | Specification                                                                                                     |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | SQLE Rabbit pAb                                                                                                   |
| Remarks/Alias    | SQLE                                                                                                              |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 6713                                                                                                              |
| GeneID           | 6713                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human SQLE (NP_003120.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q14534                                                                                                            |
| Protein Weight   | 64kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |