

Anti-Sqrdl Polyclonal Antibody For Wb, Elisa - Q9Y6N5

Item Number: CM0026197



Introduction

Rabbit polyclonal antibody against human SQRDL (sulfide:quinone oxidoreductase), validated for WB and ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial sulfide oxidation research.

[Learn More](#)

Parameter	Value
Product Name	SQRDL Rabbit pAb
Remarks/Alias	SQR; SQRDL; CGI-44; PRO1975
Species	Human
Gene ID	58472
Immunogen	Recombinant fusion protein corresponding to amino acids 211-450 of human SQRDL (NP_067022.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6N5
Protein Weight	50kDa
Shipping	Ice bag