


Anti-Sqs/Fdft1 Rabbit Monoclonal Antibody For Wb, Elisa - P37268

Item Number: CM0010313



Introduction

Rabbit monoclonal antibody targeting human SQS/FDFT1 (squalene synthase). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for cholesterol metabolism and lipid biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	SQS/FDFT1 Rabbit mAb
Remarks/Alias	SS; SQS; DGPT; ERG9; SQSD; SQS/FDFT1
Species	Human
GeneID	2222
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human SQS/FDFT1 (P37268)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37268
Protein Weight	48kDa
Shipping	Ice bag