


Anti-Sqs/Fdft1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P37268

Item Number: CM0025925



Introduction

Rabbit polyclonal antibody against human SQS/FDFT1 (Squalene synthase), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for cholesterol biosynthesis studies. Immunogen: recombinant protein (aa 1-260).

[Learn More](#)

Parameter	Value
Product Name	SQS/FDFT1 Rabbit pAb
Remarks/Alias	SS; SQS; DGPT; ERG9; SQSD; SQS/FDFT1
Species	Human
GeneID (Human)	2222
GeneID	2222
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human SQS/FDFT1 (NP_004453.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37268
Protein Weight	48kDa
Shipping	Ice bag