


Anti-Sqstm1/P62 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13501

Item Number: CM0025834



Introduction

Anti-SQSTM1/p62 polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat SQSTM1/p62 (48kDa). Key tool for autophagy, ubiquitin signaling, and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	SQSTM1/p62 Rabbit pAb
Remarks/Alias	p60; p62; A170; DMRV; OSIL; PDB3; ZIP3; p62B; NADGP; FTDALS3; SQSTM1/p62
Species	Human
GeneID (Human)	8878
GeneID	8878
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-440 of human SQSTM1/p62 (NP_003891.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13501
Protein Weight	48kDa
Shipping	Ice bag