


Anti-Sqstm1/P62 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13501

Item Number: CM0016923



Introduction

Rabbit polyclonal antibody against human SQSTM1/p62, validated for WB, IF/ICC, IP, ELISA with cross-reactivity to mouse and rat. Ideal for autophagy, ubiquitin signaling, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	SQSTM1/p62 Rabbit pAb
Remarks/Alias	p60; p62; A170; DMRV; OSIL; PDB3; ZIP3; p62B; NADGP; FTDALS3; SQSTM1/p62
Species	Human
GeneID (Human)	8878
GeneID	8878
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13501
Protein Weight	48kDa
Shipping	Ice bag