


Anti-Sr-B2/LimpII Polyclonal Antibody For Wb, If/Icc, Elisa - Q14108

Item Number: CM0016245



Introduction

Rabbit polyclonal antibody targeting human SR-B2/LIMP2 (LIMP2/LGP85/CD36L2). Validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Suitable for lysosomal receptor and enterovirus 71 entry studies.

[Learn More](#)

Parameter	Value
Product Name	SR-B2/LIMP2 Rabbit pAb
Remarks/Alias	AMRF, EPM4, LGP85, CD36L2, HLG85, LIMP-2, LIMP2, SR-B2, SR-B2/LIMP2
Species	Human
GeneID (Human)	950
GeneID	950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human SR-B2/LIMP2 (NP_005497.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14108
Protein Weight	54kDa
Shipping	Ice bag