


Anti-Sr-B2/LimpII Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14108

Item Number: CM0010008



Introduction

Rabbit monoclonal antibody targeting human SR-B2/LIMP2 (lysosome membrane protein 2), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for studying lysosomal GBA1 targeting and enterovirus 71 entry.

[Learn More](#)

Parameter	Value
Product Name	SR-B2/LIMP2 Rabbit mAb
Remarks/Alias	AMRF; EPM4; LGP85; CD36L2; HLG85; LIMP-2; LIMP2; SR-BII; SR-B2/LIMP2
Species	Human
Gene ID (Human)	950
Gene ID	950
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SR-B2/SR-B2/LIMP2 (Q14108).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14108
Protein Weight	54kDa
Shipping	Ice bag