

Anti-Sr-Bi Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wtv0

Item Number: CM0017764

Introduction

Rabbit polyclonal antibody against human SR-BI (CLA1, CD36L1), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for investigating cholesterol homeostasis, HCV entry, and scavenger receptor biology.

[Learn More](#)

Parameter	Value
Product Name	SR-BI Rabbit pAb
Remarks/Alias	CLA1; SRB1; CLA-1; SR-BI; CD36L1; HDLCQ6; HDLQTL6
Species	Human
GeneID (Human)	949
GeneID	949
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SR-BI (NP_005496.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTV0
Protein Weight	61kDa
Shipping	Ice bag