

## Anti-Sr Protein Repeat Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q13247

Item Number: CM0012357



### Introduction

This rabbit monoclonal antibody recognizes SR protein repeat (SRSF6/SRP55) in human, mouse, and rat. Validated for WB, IHC-P, IP, and ELISA. Molecular weight ~40 kDa. Suitable for splicing research and disease model studies.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | SR protein repeat Rabbit mAb                                                                                   |
| Remarks/Alias    | B52; SFRS6; SRP55; HEL-S-91; SR protein repeat                                                                 |
| Species          | Human                                                                                                          |
| GeneID (Human)   | 6431                                                                                                           |
| GeneID           | 6431                                                                                                           |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 100-200 of human SR protein repeat (Q13247) |
| Source           | Rabbit                                                                                                         |
| Category         | Monoclonal Antibodies                                                                                          |
| Application      | WB, IHC-P, IP, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                              |
| SWISS            | Q13247                                                                                                         |
| Protein Weight   | 40kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                        |