


Anti-Sra1 Polyclonal Antibody For Wb, Elisa - Q9Hd15

Item Number: CM0021393



Introduction

Anti-SRA1 rabbit polyclonal antibody for WB and ELISA. Targets human SRA1 (SRAP), a transcriptional coactivator involved in steroid receptor regulation and tumorigenesis. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	SRA1 Rabbit pAb
Remarks/Alias	SRA; SRAP; STRAA1; pp7684; SRA1
Species	Human
Gene ID (Human)	10011
Gene ID	10011
Immunogen	Synthetic Peptide: a synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRA1 (NP_001030312.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9HD15
Protein Weight	24kDa
Shipping	Ice bag